

Instrument Security Procedures

Model:

Fluke 411D, 416D, 424D

Product Name:

Distance meter 411D, 416D, 424D

Instrument Description:

The Fluke 411D, 416D, and 424D are Laser Distance meter tools.

Memory Description:

The Fluke 411D, 416D, and 424D has the following memory devices:

1. Microcontroller ROM or Flash for system software.
2. EEPROM for calibration.

Memory Cleaning Instructions:

1. The operating code (instrument firmware) stored in 411D, 416D, and 424D
2. Calibration constants stored in 411D, 416D, and 424D can be read using special remote interface commands