

Î

## **Recommended pharmaceutical solutions**

#### **Power quality loggers and analyzers**

#### Best suited for:

- Troubleshooting and preventing problems in power generation, transmission and distribution systems
- Conducting load studies to reduce energy usage • Use with motors, pumps, generators, transformers, compressors and electrical panels



#### **Field metrology wells**

#### **Best suited for:**

- Automated comparison calibrations of RTDs and thermocouples (-95°C to 660 °C)
- Transmitter loop calibrations
- Thermal switch tests
- 21 CFR Part 11 compliance with TQSoft Thermal Validation Software and the Fluke Calibration 1586A Super DAQ Precision Temperature Scanner
- Pairs with Fluke 754 Documenting Process Calibrator



#### **Multifunction process calibrators**

#### Best suited for:

- Calibrating and maintaining process instrumentation and control systems
- Troubleshooting 4–20 mA loops, temperature,
- pressure, frequency and transmitters in one tool Performing HART transmitter adjustment



### **Thermal imagers Best suited for:**

#### 0 HVAC inspections

- Identifying "hot spots" in motors, pumps and compressors
- Detection of gas leaks in refrigeration systems







Fluke Connect

• Simplifying preventive

data into software

maintenance programs

• Seamlessly integrating tool

• Diagnosing asset health in

• Cloud-based data sharing

Best suited for:

the field

#### **Indoor air quality**

**Best suited for:** • Troubleshooting and monitoring indoor air quality • Optimizing HVAC ventilation settings





# **Pharmaceutical Solutions**

Pharmaceutical manufacturing requires premium quality tools that provide accurate measurements and keep critical processes up and running. Whether you're maintaining plant uptime, documenting product quality, or calibrating measurement devices, you need accurate tools that you can count on for the bench and the manufacturing floor.

With solutions from Fluke for electrical, multifunction, mA loop, temperature, pressure and flow applications, you can find the right tools for the challenges you face every day.

Fluke. Keeping your world up and running.

**Fluke Corporation** PO Box 9090, Everett, WA 98206 U.S.A. Fluke Europe B.V. PO Box 1186, 5602 BD Eindhoven, The Netherlands

For more information call: In the U.S.A. (800) 443-5853 or Fax (425) 446-5116 In Europe/M-East/Africa +31 (0) 40 2675 200 or Fax +31 (0) 40 2675 222 In Canada (800)-36-FLUKE or Fax (905) 890-6866 From other countries +1 (425) 446-5500 or Fax +1 (425) 446-5116 Web access: http://www.fluke.com

©2015, 2017 Fluke Corporation. Specifications subject to change without notice. Printed in U.S.A. 8/2017 4264583c-en

**Best suited for:** • Pressure calibration of transmitters, gauges, switches and custody transfer meters Maintenance of almost any pressure device using

Pressure

calibrators

- pressure sourcing and mA measurements Calibrating P/I instruments
- with simultaneous pressure and current measurement





### **Building Control Technician**

#### **Issues faced:**

- Troubleshooting building environmental systems with immediate response
- Maintaining uptime and production area cleanliness
- Indoor air quality • Safety systems testing

Fluke solutions:

Fluke Connect

▲ Loop calibration tools

• Thermal imagers

### **Quality Manager**

- **Issues faced:**
- Complying with quality regulations
- Maintaining documentation of processes
- Fluke Connect

Fluke solutions:

- A Multiprocess calibration tools
- Loop calibration tools
- A Temperature calibration tools
- Pressure calibration tools
- + Indoor air quality
- Thermal imagers

### **Reliability Engineer**

#### Issues faced:

- Specifying equipment to meet regulations and requirements
- Preventing downtime

### Fluke solutions:

- Fluke Connect A Multiprocess calibration tools
- Loop calibration tools
- A Temperature calibration tools
- Pressure calibration tools
- Thermal imagers

## Fluke pharmaceutical solutions

Learn how Fluke solutions can save you time and money while increasing productivity at www.fluke.com.

Multiprocess Calibration Tools	Process Calibration Loop Tools		Process Calibration Temperature Tools				Calibration	Pressure	e Tools	Electrical Test Tools		Infrared	Air Quality	Software
Documenting Multifun process calibrators process cal	Process test tools	Field metrology wells		Reference thermometers	Temperature calibrators	Pressure gauges	Pressure calibrators		Power quality analyzers	Digital multimeters	Thermal imagers professional	Indoor air quality	Maintenance software	
	773	P190A	9142	1523	724	700G	3130	729	718-16	430 Series II 17		) • • • • • • • • • • • • • • • • • • •	985	Fluke Connect*





### **Field Calibration and Instrumentation Technician**

#### **Issues faced:**

- Preventing downtime Troubleshooting process
- interruptions
- Testing safety systems
- Maintaining product
- accuracy for regulations
- and quality

### Fluke solutions:

- Fluke Connect
- ▲ Loop calibration tools
- Temperature calibration tools
- Pressure calibration tools
- + Indoor air quality
- Electrical test tools
- Thermal imagers

### Validation Technician

#### **Issues faced:**

- Thermal validations
- and studies
- Authoring and executing
- validation protocols, and
- writing reports of validation results

#### Fluke solutions:

Fluke Connect

- A Multiprocess calibration tools
- A Temperature calibration tools

### Lab Calibration

#### **Issues faced:**

- Making accurate
- measurements
- Documenting procedures
- Increasing throughput

#### Fluke solutions:

- A Multiprocess calibration tools
- Loop calibration tools
- Temperature calibration tools
- Pressure calibration tools