

Instrument Security Procedures

Model: 709 and 709H

Product Name: 709 Precision Loop Calibrator, 709H Precision Hart Loop Calibrator

Instrument Description: 4-20 mA loop calibrator. The H version has HART communication capability.

Memory Description:

Microprocessor:

128 Kbytes flash

4 Kbytes EEPROM

Microprocessor RAM:

4 Kbytes

Data logging EEPROM:

256 K (709H only)

Memory Cleaning Instructions:

Use Programming tool to erase microprocessor flash and EEPROM.

Use Command from keypad or PC to erase data logging memory (709H only)