

FLUKE EX DEMO CLEARANCE

FLUKE.

Fluke is offering thousands of dollars in savings on selected Fluke ex demo clearance tools. Limited stock is available. Items sell out fast so get in QUICK!

Prices available until stock lasts.
Prices listed are Recommended Trade Price excluding GST

Deals only available via Fluke Australia direct. Purchase your Fluke product by filling out the order form below and sending via fax 61 2 8850 3300 or email igsales@fluke.com.au

<p>Fluke Ti32 60Hz Infrared Camera + Wide1 or Tele1 Lens + Visor + Tripod</p>  <p>Was \$14,890 Now \$9,790</p> <p>SOLD OUT</p> <p>1 LEFT</p>	<p>Fluke Ti125 9Hz Infrared Camera</p>  <p>Was \$8,040</p> <p>SOLD OUT</p> <p>1 LEFT</p>	<p>Fluke 1555 10 kV Insulation Resistance Tester</p>  <p>SOLD OUT</p> <p>3 LEFT</p>	<p>Fluke T5-1000 Electrical Tester</p>  <p>Was \$277 Now \$215</p> <p>3 LEFT</p>	<p>Fluke 62 Max+ Infrared Thermometer</p>  <p>Was \$267 Now \$145</p> <p>5 LEFT</p>
<p>Fluke 1630 Earth Ground Clamp Meter</p>  <p>Was \$2,600 Now \$1,465</p> <p>SOLD OUT</p> <p>1 LEFT</p>	<p>Fluke 375 True RMS Clamp Meter</p>  <p>Was \$755 Now \$425</p> <p>SOLD OUT</p> <p>1 LEFT</p>	<p>Fluke 376 True RMS Clamp Meter</p>  <p>Was \$885 Now \$650</p> <p>2 LEFT</p>	<p>Fluke 381 Remote Display True RMS Clamp Meter</p>  <p>Was \$950 Now \$700</p> <p>1 LEFT</p>	<p>Fluke 902 True RMS HVAC Clamp Meter</p>  <p>Was \$595 Now \$325</p> <p>2 LEFT</p>
<p>Fluke 709H Precision Loop Calibrator with HART + Fluke 709H/Track HART Communications Software</p>  <p>Was \$2,481 Now \$1,895</p> <p>SOLD OUT</p> <p>1 LEFT</p>	<p>Fluke 80PK-8 Pipe Clamp Temperature Probe</p>  <p>Was \$335 Now \$170</p> <p>6 LEFT</p>	<p>Fluke 51 II Single Input Digital Thermometer</p>  <p>Was \$539 Now \$375</p> <p>1 LEFT</p>	<p>Fluke CNX C3000 Premium Modular Tool Bag</p>  <p>Was \$175 Now \$75</p> <p>1 LEFT</p>	

Ex Demo Clearance Order Form

FLUKE.

SOLD OUT
Fluke Ti220 Infrared
Camera or Tele1 Lens
+ Tripod
\$9,790.00 exc. GST

SOLD OUT
Fluke In...
\$4,550.00 exc. GST

SOLD OUT
Fluke 1555 V Insulation
Resistor
\$4,000.00 exc. GST

Fluke T5-1000
Electrical Tester
\$215.00 exc. GST

Fluke 62 Max
Infrared Thermometer
\$145.00 exc. GST

SOLD OUT
Fluke 375
True RMS Meter
\$320.00 exc. GST

Fluke 376
True RMS Clamp Meter
\$650.00 exc. GST

Fluke 381 Remote Display
True RMS Clamp Meter
\$700.00 exc. GST

Fluke 902
True RMS HVAC Clamp Meter
\$325.00 exc. GST

SOLD OUT
Fluke 1620
Ear... Meter
\$1,400.00 exc. GST

SOLD OUT
Fluke 700... on Loop
Calibration HART
Communications Software
\$1,895.00 exc. GST

Fluke 80PK-8 Pipe Clamp
Temperature Probe
\$170.00 exc. GST

Fluke 51 II Single Input
Digital Thermometer
\$375.00 exc. GST

Fluke CNX C3000 Premium
Modular Tool Bag
\$75.00 exc. GST

CUSTOMER DETAILS

Please complete your full details below.

Name:	
Job Title:	
Company:	
Address:	
Suburb:	State:
Postal Code:	Telephone:
Fax:	Email:

PAYMENT DETAILS

CREDIT CARD: VISA MASTER CARD *AMEX not accepted

Name on card: _____

Card No.: _____

Expiry Date: _____

Signature: _____

Bill my Credit Card for the total amount of: _____ + GST

CVC No. (Last 3 digits on back of card or 4 digits on front): _____

Please fax/email your completed order form to 02 8850 3300 or igsales@fluke.com.au

TERMS AND CONDITIONS

Sponsored by Fluke Australia. Fluke Australia Pty Ltd conditions apply. Prices are in AUD excluding GST. Savings percentages have been rounded to the nearest whole number. Prices good only in Australia, New Zealand and Pacific Island regions through Fluke Australia direct with Credit Card payment (excluding AMEX). Promotion cannot be combined with any other Fluke offer. Items will be sold in a "as is" status. Products sold on a first in first served basis. Items will come with new 6 month warranty. Promotional prices are valid until stocks last. Please allow 4-6 weeks for delivery.

ORDER DETAILS

To confirm your order please ensure the above details are complete in full and you fax or scan this form and send to **Fluke Australia on FAX 02 8850 3300 or EMAIL igsales@fluke.com.au**