

## Instrument Security Procedures

**Model:**

CNX 3000

**Product Name:**

CNX 3000 Wireless Multimeter

**Instrument Description:**

CNX 3000 Wireless Multimeter

**Memory Description:**

Fluke CNX 3000 has one memory device:

- U6, micro controller, TI MSP430F5435A, 16-BIT, 1.8 V-3.6V, 25 MHZ, 192 KB FLASH, 16 KB RAM

The user-controlled storage operation (other than calibration constants) is the storage of Min Max and Average.

**Memory Cleaning Instructions:**

The memory is cleared when the instrument is powered off.